BRAVO TD III

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Instruction for Use

Mode d'emploi

Bedieningshandleiding

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Manual de instrucciones

Instruções

BRAVO TD III Labormotor

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BRAVO TD III Laboratory Motor


Bravo TD III is used for minor polishing and grinding works with hand- and/or angle-pieces in the practice and laboratory field and **is intended for extraoral use only**.

Initiation

Please make sure that the Bravo TD II device is standing on a fixed, even surface. The connectors on the back are:

- (1) power cable
- (2) foot control
- (3) micro-fuse, 1.6 A

Plug in the power cable and the foot control and connect the device to the power socket.


At the front of the device you find:

- (a) on-/off-switch
- (b) revolution speed control
- (c) reverser clockwise-/counter-clockwise rotation
- (d) motor connection

Connect the motor to the device. Insert plug and turn clockwise in a 45° angle, mount hand-and angle-piece to the motor, switch-on the device, select the desired rotation speed and activate foot control. The device is equipped with an overload protection. The motor stops in case of overload. For re-starting, the foot control has to be actuated anew. Keep objects, like rags, cloths and hair away from the grinding head when in operation.

In order to guarantee a long motor life, please don't insert freshly oiled hand- and angle-pieces. Allow the excess to drip off. For this, put hand- and angle-pieces upright on an absorbent pad. Don't change the direction of rotation during operation. This may only be done when the motor is out of operation.

Start motor only with attached hand- or angle-piece.

Warranty

The maker guarantees that the contractual product is free of defects and features all warranted and/or agreed characteristics. The warranty period begins with the day of transfer of risks or delivery of the device and lasts for twelve months. Excluded from warranty are: defects and/or damages due to e.g. operational wear and tear, normal wearout, improper use, operation faults and negligence of the customer when handling the device, connection to inappropriate power sources, operation with wrong power voltage, force majeure, such as fire, lightning, humidity, etc. In case of justifiable defects, the maker has the right to repair or compensation delivery or reduction of the purchase price. The customer may only claim cancellation of the contract after two unsuccessful repairs by the maker.

Technical Data

Revolution speed	1,000 – 30, rpm
Motor capacity	max. 40 W/three-phase
Motor dimensions	dia. 22 mm x 114 mm length
Power supply casing	80mm x 200 mm x 140 mm
Weight	1.2 kg
Supply voltage	230 V (50 Hz)
Art.-no.	452 106

The set consists of:
control box, micro motor, foot control and power cable.

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HÄGER & WERKEN GmbH & Co. KG